



Kineticopro[®]
Total Water Care

Kineticopro Water Test Kit Reference Guide

1. You will receive a box labeled 'Water Sample Test Kit' from Kineticopro. This will include a return shipment label once you have completed the instructions inside.

2. Inside the box you will find 2 bottles. Please remove these items from the box.

NOTE: We have discontinued the ph strips

3. If you are a Service Provider and have received the "Water Sample Collection." sheet, please fill it out.

4. The 1st bottle, labeled '**Unfiltered Water**', will be used to capture a sample of water from water source that IS NOT currently being treated. (ex: bathroom sink, kitchen sink, etc.).

5. The 2nd bottle, labeled '**Filtered Water**', will be used to capture a water sample from a treated water source, if applicable (ex: water from fountain beverage machine, filtered drink water station, etc.).

6. When filling up the water bottles please be sure to follow the instructions provided on the side of each water bottle, which can also be found inside the box.

7. After **water bottles are filled all the way up**, please ensure lids are screwed on tightly and put back into the box provided. Enclose (or email) the "Water Sample Collection" sheet. Then adhere return label on the outside of the box and call UPS to make pick-up arrangements OR drop off at your nearest UPS store.

All bottles need to be returned to the following address:

Kineticopro Incorporated
3980 Lakefield Ct.,
Suwanee GA, 30024

Please email the Water Sample Collection sheet to:
commercial@kinetico.com

1.



2.



3.



4.



5.



6.



7.



1.800.321.5022
www.Kineticopro.com
WaterCare@Kineticopro.com